

Comparison Chart

Lateral Turning Equipment

This list provides a starting point for assessors who are looking at equipment options from a range of suppliers.

This list includes some options requested through Enable New Zealand and does not include all options available in NZ. Information is current as of the date shown in the footer and subject to change.

Product	Image	Safe working load (kg)	Supplier and Product Details
TurnAll (Levabo)		Standard: 225 Bariatric: 300	Active Healthcare Details: - To be used between the bed base and mattress 30 degrees of lateral tilt with 5 degrees opposite side support - Cycle time: 30min, 60min and 90min intervals - Side isolation: Can isolate turning to just one side via an adapter (separate item in quote) Dimensions: - W850mm x L2000mm - Bariatric mattress: W1100mm x L2000mm Supplier link: https://activehealthcare.co.nz/products/pressure-care/pressure-relief-aids/turnall/

Toto Lateral Turning System		250	Cubro Details: • To be used between the bed base and mattress • 25 degrees of lateral tilt • Cycle time: 30min, 60min, 120min, 180min, 240min • Side isolation: turning can be isolated to a single side if required • Recommended to be used with side rails • Cubro indicates it can be used for single bed (900mm wide) and if the unit and client are positioned well (in the centre of the bed and mattress), it can be used on a super single bed (1100mm wide) Dimensions: • W750mm x L1800mm Website link: • https://www.cubro.co.nz/healthcare-products/pressure-care/lateral-turning-systems/toto-lateral-system • Cubro have videos on their website discussing/demonstrating the Toto turning system
--	---	------------	---

Vendlet V5S Turning System		V5S: 200 Bariatric V5S+: 400 Minimum user weight: 40kg	Cubro Details: • To be used on top of a mattress • Electronic turning sheet • Adjustable turning speeds: 50%, 75% or 100% • Electric height adjustable rails • Bed width is not critical as extra wide turning sheets available – check with Cubro (Single, super single, double) • Check with Cubro whether the Vendlet is compatible with current bed Dimensions: • W75mm x L2000mm on either side of the bed Website Link: • https://www.cubro.co.nz/healthcare-products/moving-and-handling/liftingand-transfer-aids/vendlet-v5s-turning-system
TurnAid Turning System		350 Minimum user weight: 40kg	Cubro Details: • To be used on top of a mattress • Electronic turning sheet • Adjustable turning speeds: 50%, 75% or 100% • Sides fold out and down (manual action) • Can be used on a single or super single sized bed • Foot and/or hand control • Compatible with extra low beds • Check with Cubro whether the TurnAid is compatible with current bed dimensions Dimensions: • Single: W900mm x L1900mm x H570mm

			• Super Single: W1100mm x L1900mm x H570mm Website Link: • https://www.cubro.co.nz/healthcare-products/moving-and-handling/liftingand-transfer-aids/turnaid-safe-turning-system • Cubro have videos on the TurnAid webpage discussing/demonstrating the TurnAid system
--	--	--	---