



**Enable**  
NEW ZEALAND®

## **Complex Lateral Turning Equipment**

Clinical Best Practice to Optimise Pressure  
Care and Safety

Presented by:

Dheshini Naidoo, Jolene Young, Pauline Lazarus and Lisa Brown  
Clinical Service Advisors, Enable New Zealand

# Karakia timatanga

Kia hora te marino

Kia whakapapa pounamu te moana

Hei huarahi mā tatou i te rangi nei

Aroha atu, aroha mai

Tātou i a tatou katoa

Hui e, Tāiki e!

May peace be widespread

May the sea be like greenstone

A pathway for us all this day

Let us show respect for each other

For one another

Bind us all together!

# Introductions



**Jolene Young**

Occupational Therapist  
Wheelchairs, seating, equipment and  
housing  
Phone: (03) 740 3912  
Email: [jolene.young@enable.co.nz](mailto:jolene.young@enable.co.nz)



**Lisa Brown**

Occupational Therapist  
Housing and equipment  
Phone: (03) 740 3911  
Email: [lisa.brown@enable.co.nz](mailto:lisa.brown@enable.co.nz)



**Pauline Lazarus**

Clinical Services Advisor  
Physiotherapist  
Phone: (07) 843-9181



**Dheshini Naidoo**

Clinical Services Advisor  
Occupational Therapist  
Phone: (07) 981 1369

# Objectives

- Build confidence in prescribing complex turning equipment.
- Share current international guidelines, and findings from our literature review.
- Provide clinical insights from recent equipment requests.
- Apply the Internal Classification of Function (ICF) framework (WHO, 2001) to guide assessment and analysis .
- Discuss clinical reasoning processes to support equipment prescription.
- Present a comparison chart with available equipment options.
- Facilitate a forum for discussion and Q&A.

# Turning Systems Available in Aotearoa

| Automated                                   |                                                                                       | Requires x1 carer support |                                                                                       |
|---------------------------------------------|---------------------------------------------------------------------------------------|---------------------------|---------------------------------------------------------------------------------------|
| TurnAll (Levabo)<br>(Active Healthcare n.d) |   | TurnAid<br>(Cubro n.d)    |   |
| Toto<br>(Cubro n.d)                         |    | Vendlet<br>(Cubro n.d)    |  |
| Sapphire Comfort Turn<br>(Arjo n.d)         |  |                           |                                                                                       |

# International Guidelines

(International Guideline, 2025)

## Positioning the individual in bed – Implementation Considerations

Some full body support surfaces offer **automated lateral turning**. These are available in integrated bed/mattress systems or as devices that can be positioned on an existing full body support surface.


There is **limited evidence on the efficacy of these devices on PI occurrence**; however, they **may be appropriate** for individuals **requiring frequent repositioning** (e.g., in critical care settings) or for individuals who have **limited access to assistance to reposition** (e.g., in home care settings) (23, 88, 89).

# Our Literature Search

| Journal Article                                                                                                                                                                                | Type of research & participants                                                                                                                                                                                             | Intervention                                                                                                                                                                                                                                                                                                 | Outcome                                                                                                                                                                                                                                                                                                     | Limitations                                                                                                                                                                                                                         |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>An automatic repositioning system to prevent pressure ulcers: a case series</p> <p><i>Knibbe et al 2017</i><br/><i>British Journal of Nursing</i></p>                                       | <ul style="list-style-type: none"> <li>Case series</li> <li>Residents in 2 Residential Care Facilities.</li> </ul>                                                                                                          | <ul style="list-style-type: none"> <li>4 weeks use of <b>Vendlet</b> V5.</li> <li>Re-positioning protocol was individualised for each person.</li> </ul> <p>Measurements –</p> <ul style="list-style-type: none"> <li>Skin checks 2x/week.</li> <li>Self-report of carer musculoskeletal symptoms</li> </ul> | <ul style="list-style-type: none"> <li><b>Only 1 nurse</b> was required to reposition using the Vendlet.</li> <li><b>Fewer carer musculoskeletal symptoms</b> were reported.</li> <li>The authors indicated that no new PIs developed and that 4 pre-existing PI's healed during the experiment.</li> </ul> | <ul style="list-style-type: none"> <li><b>n = 13</b> residents</li> <li><b>Potential for bias</b> as no blinding, randomisation and control group.</li> <li>Outcome measures not proven to have reliability or validity.</li> </ul> |
| <p>Real-world experience of using the <b>ToTo</b> Lateral Turning System in a busy spinal unit: benefits and tips in practice.</p> <p><i>Rebecca Warren and Gemma Kerr 2019, Wounds UK</i></p> | <ul style="list-style-type: none"> <li>Product evaluation at Rehab centre (15 acute and 31 rehab beds).</li> <li>Used with patients unable to independently turn.</li> <li>Not used if there was an existing PI.</li> </ul> | <ul style="list-style-type: none"> <li>Use of <b>ToTo</b> lateral turning system as standard practice for reducing risk of pressure damage.</li> <li>Used for suitable and willing patients.</li> </ul>                                                                                                      | <ul style="list-style-type: none"> <li>Overall <b>positive feedback</b> by patients and staff.</li> <li>ToTo much <b>less disruptive on sleep</b> vs manual repositioning.</li> <li>Sometimes improved <b>comfort</b>.</li> <li><b>Reduces staff time</b> taken to reposition people.</li> </ul>            | <ul style="list-style-type: none"> <li>Article is a narrative/<b>opinion</b>.</li> <li>The authors indicate the ToTo turning system is clinically effective, however they have <b>not referenced</b> this statement.</li> </ul>     |

# Review of Solutions Proposed

## Request Number



## Diagnoses & situation

### TurnAll & Toto:

- Recent C3/4 Spinal Cord Injury (SCI)
- SCI years ago, on bed rest due to a pressure injury.
- DSS Funding – predominantly neurological conditions and 50% had an active pressure injury.


### Vendlet & TurnAid:

- Clients fully dependent for cares, bed mobility & transfers
- Diagnoses usually Tetraplegia, Brain Injury or other neurological condition.

# Goals When Prescribing Equipment

| TurnAll/Toto                                              | Vendlet/TurnAid                                |
|-----------------------------------------------------------|------------------------------------------------|
| Night pressure redistribution/offload PIs(15)             | Reduce demand on carers (9)                    |
| Reduce need for manual handling and carers overnight (14) | More comfortable than manual repositioning (5) |
| Reduce sleep disturbance & improve sleep (4)              | Meet repositioning needs (2)                   |
| Avoid shear (2)                                           | Repositioning with just 1 carer (1)            |
|                                                           | Reduce shear (1)                               |

# Outcomes of trial – TurnAll & Toto




## Successful

- No longer waking for repositioning
- Improved sleep
- Pressure injuries improved
- effective for individuals with low body weight
- effective for individuals with shoulder pain affecting bed mobility.

## Unsuccessful

- Body weight over 120kg – felt unsafe and unable to use the whole mattress
- Did not tolerate the mattress movement
- Did not want to sleep in supine.

# Outcomes of trials – Vendlet & TurnAid



## Successful

- Turnaid more robust than Vendlet
- Reduced carer reach using Vendlet
- Vendlet – can lower rail at one end and leave it up at the other end
- worked for people with higher body weight.


## Unsuccessful

- Vendlet motor hits floor if used on ultra low beds.

# Our Review – Summary

- Lateral turning devices could also be considered when pressure injuries are not resolving.
- The TurnAll and Toto may also improve the persons sleep, assist with their nighttime comfort, and allow their carers to sleep through the night.
- The TurnAid and Vendlet may reduce overnight carers from 2 to 1 person, reduce carer musculoskeletal strain and improve the persons comfort when being re-positioned.
- Ongoing skin checks remain essential if equipment is put in place long term

# Assessment using the ICF framework



(WHO, 2001)

# Clinical Reasoning

- What are the client and whānau/carer goals?
- What am I aiming to achieve or improve for the client?
- What are the clients current challenges/barriers?
- What supports and equipment have been considered and discounted?
  - supports within the home
  - low cost equipment
  - is the current mattress meeting their need.
- Proposed equipment – How will the complex lateral turning equipment benefit the client, and what positive changes will it bring?
- Does the solution meet an essential need and represent value for money?

# Comparison Chart

| Product name                | Image                                                                               | Safe working load (kg)                 | Supplier and Product Details                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-----------------------------|-------------------------------------------------------------------------------------|----------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| TurnAll (Levabo)            |    | 225kg<br><br>Bariatric mattress: 300kg | <p>Supplier: Active Healthcare</p> <p>Product Details:</p> <ul style="list-style-type: none"> <li>To be used between the bed base and mattress</li> <li>30 degrees of lateral tilt with 5 degrees opposite side support</li> <li>Cycle time: 30min, 60min and 90min intervals</li> <li>Side isolation: Can isolate turning to just one side via an adapter (separate item in quote)</li> </ul> <p>Dimensions:</p> <ul style="list-style-type: none"> <li>W850mm x L2000mm</li> <li>Bariatric mattress: W1100mm x L2000mm</li> </ul> <p>Website link: <a href="https://activehealthcare.co.nz/products/pressure-care/pressure-relief-aids/turnall/">https://activehealthcare.co.nz/products/pressure-care/pressure-relief-aids/turnall/</a></p> <p>Active Healthcare video link: <a href="https://activehealthcare.co.nz/webinar/">https://activehealthcare.co.nz/webinar/</a> (password: active)</p> |
| Toto Lateral Turning System |  | 250kg                                  | <p>Supplier: Cubro</p> <p>Product Details:</p> <ul style="list-style-type: none"> <li>To be used between the bed base and mattress</li> <li>25 degrees of lateral tilt</li> <li>Cycle time: 30min, 60min, 120min, 180min, 240min</li> <li>Side isolation: turning can be isolated to a single side if required</li> <li>Recommended to be used with side rails</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |

# Future Research and Clinical Questions

The Guideline Governance Group noted the following future clinical questions for potential/future address in the guideline:

- *Should an automated repositioning device (e.g. lateral tilt device) versus manual repositioning be used to prevent pressure injury occurrence in individuals at risk?*
- *Should repositioning turn systems versus not using a turn system be used to prevent pressure injury occurrence in individuals at risk?*

(International Guideline, n.d.)

# Key Points

- The International Guidelines state limited evidence of the efficacy of automated lateral turning device on pressure injury occurrence and highlights the need for further research in this area.
- Our review of complex turning equipment demonstrated effectiveness in pressure injury management, reducing carer burden, and improving comfort and sleep.
- A comprehensive assessment and sound clinical reasoning are essential.
- Refer to the new 2025 PI clinical practice guideline for latest evidence-based practice.

# Feedback

Help us to support you with education opportunities you want and need

**These presentations are to support your learning and development.**



We want to hear what areas you would like to explore.

We welcome any comments or feedback about today's session.



<https://www.surveymonkey.com/r/QP56PBC>

Scan the QR code or click the link to answer

# Questions / Contact details



## **Email to gain clinical input: –**

ACC Assessors: [acc.advisor@enable.co.nz](mailto:acc.advisor@enable.co.nz)

DSS Assessors: [CSA@enable.co.nz](mailto:CSA@enable.co.nz)

## **Email for general queries –**

[enable@enable.co.nz](mailto:enable@enable.co.nz)

**Phone – 0800 362 253**

# Closing Karakia

**Kia whakairia te tapu  
Kia wātea ai te ara  
Kia turuki whakataha ai  
Kia turuki whakataha ai  
Haumi e, Hui e, Tāiki e!**

**Restrictions are moved aside  
So the pathways are clear  
To return to everyday activities**

# References

- Active Healthcare. (n.d.). TurnAll pressure relief aid. <https://activehealthcare.co.nz/products/pressure-care/pressure-relief-aids/turnall/>
- Arjo. (n.d.). Sapphire Comfort Turn. <https://www.arjo.com/en-gb/products/pressure-ulcer-prevention---pup/alternating-pressure/sapphire-comfort-turn/>
- Arjo. (n.d.). Pressure injury/ulcer prevention: Understanding the science of support surfaces [Webinar]. Arjo. <https://www.arjo.com/en-au/knowledge/global-webinars/pressure-injuryulcer-prevention-understanding-the-science-of-support-surfaces/>
- Cubro. (n.d.). Cubro healthcare solutions. <https://www.cubro.co.nz/>
- International Guideline. (2025, February 25). Repositioning for pressure injury prevention. [https://internationalguideline.com/repositioning#s2\\_info](https://internationalguideline.com/repositioning#s2_info)
- International Guideline. (n.d.). Future research. <https://internationalguideline.com/future-research>
- Knibbe, H., et al. (2017). An automatic repositioning system to prevent pressure ulcers: A case series. British Journal of Nursing, 26(15), S6–S12.
- Warren, R., & Kerr, G. (2019, October 31). Real-world experience of using the ToTo Lateral Turning System in a busy spinal unit: Benefits and tips in practice. Wounds UK, 15(5). <https://wounds-uk.com/journal-articles/real-world-experience-of-using-the-toto-lateral-turning-system-in-a-busy-spinal-unit-benefits-and-tips-in-practice/>
- World Health Organization. (n.d.). International classification of functioning, disability and health (ICF). <https://www.who.int/standards/classifications/international-classification-of-functioning-disability-and-health>